



# Simplified Automation of the GMP production of $^{68}\text{Ga}$ -labelled peptides

David Goblet, Cristiana Gameiro, Neva Lazarova - IBA SA, Louvain-La-Neuve/BE


## Introduction

Optimized automated process for the GMP production of  $^{68}\text{Ga}$ -DOTA-NOC and  $^{68}\text{Ga}$ -PSMA-11 using a cassette-based synthesizer (Synthera<sup>®</sup>, IBA, Louvain-la-Neuve, Belgium) in combination with a recently commercialized  $^{68}\text{Ge}/^{68}\text{Ga}$  generator (Galli Eo<sup>™</sup>, IRE-Elit, Fleurus, Belgium) have been developed in this work.

## Materials & Methods



Set-up in a standard hot cell:  
Synthera<sup>®</sup> synthesizer combined  
with IRE Galli Eo<sup>™</sup>  $^{68}\text{Ge}/^{68}\text{Ga}$   
generator



Schematic diagram depicting  
the automated synthesis of  
 $^{68}\text{Ga}$ -labelled peptides

### $^{68}\text{Ge}/^{68}\text{Ga}$ Generator elution

- 1,1 ml – 0.1 M HCl
- No pre-purification required
- Eluate sent directly to the IBA Synthera<sup>®</sup> reaction vessel

### DOTA-NOC labelling

- 50 µg peptide
- acetate buffer: 1 mL 250 mM (pH 5)
- heated at 120°C for 5 min

### PSMA-11 labelling

- 10 µg peptide
- acetate buffer: 1 mL 1.5 M (pH 4.5)
- heated at 95°C for 5 min

### $^{68}\text{Ga}$ -DOTA-NOC purification

- HLB cartridge / washed with 10 mL of water

### $^{68}\text{Ga}$ -PSMA-11 purification

- Sep-Pak<sup>®</sup> Light C18 cartridge/- washed with 10 mL of water

### $^{68}\text{Ga}$ -DOTA-NOC formulation & sterile filtration


- Elution: 1 mL of EtOH/water 65:35 v/v to pre-loaded vial (8mL saline solution) through 0.22µm filter (Millex<sup>®</sup>-GV)

### $^{68}\text{Ga}$ -PSMA-11 formulation & sterile filtration

- Elution: 2 mL of EtOH/water 1:1 v/v mix + 2 mL PBS to pre-loaded vial (6 ml PBS) through 0.22µm filter (Millex<sup>®</sup>-GV)



## Results

$^{68}\text{Ga}$ -DOTA-NOC is produced in <20 min with  $81.5 \pm 5.2$  % radiochemical yield (RCY) (decay-corrected-d.c.) and  $^{68}\text{Ga}$ -PSMA-11 is produced in 13 min with  $97.4 \pm 2.5$  % RCY (d.c.). Reported process times include generator elution and formulation. In both cases, final products show high radiochemical purity (TLC > 97 % and 99 % respectively).



### Radio-HPLC analysis of $^{68}\text{Ga}$ -DOTA-NOC (Rt: 10 min)

HPLC methods are performed with Agilent Eclipse XDB-C18 (5 µm, 4.6 x 150 mm) (reverse phase column at flow rate of 2 mL/min. Eluent: A=Water/TFA 99.9:0.1 v/v ; B=ACN/TFA 99.9:0.1 v/v Gradient: 0 min 25%B + 75%A; 4 min: 25%B + 75%A; 10 min. 32% B + 68% A; 14 min. 50%B + 50%A; 15 min. 25%B+ 75% A, 20 min. 25% B+ 75% A.



### Radio-TLC analysis of $^{68}\text{Ga}$ -DOTA-NOC (left) and $^{68}\text{Ga}$ -PSMA-11 (right)

Radio-TLC were performed with Agilent iTLC-SG silica gel plate. Eluent: Ammonium acetate 1 M and Methanol (50:50)

## Discussion & Conclusion

Automated processes for the production of both  $^{68}\text{Ga}$ -DOTA-NOC and  $^{68}\text{Ga}$ -PSMA-11 have been successfully achieved using a commercial synthesizer and a  $^{68}\text{Ge}/^{68}\text{Ga}$  generator. The labeling procedures are straightforward and efficient, thanks to the low elution volume and high purity of the generator eluate (no need for fractionation or post-elution purification).