

Cristiana Gameiro<sup>1</sup>, Ian Oxley<sup>1</sup>, Vasko Kramer<sup>2</sup>, Antero Abrunhosa<sup>3</sup>, Maria Vosjan<sup>4</sup>, Arnold Spaans<sup>4</sup>

<sup>1</sup> IBA SA, Louvain-La-Neuve/BE, <sup>2</sup> PositronPharma, Santiago/CL, <sup>3</sup> ICNAS, Coimbra/PT; <sup>4</sup> BV Cyclotron VU, Amsterdam/NL

## Introduction


Mono-product facilities are turning to multi-product to cope with decreasing prices of [<sup>18</sup>F]-FDG (FDG) and to support R&D programs. The aim is to demonstrate that in the same facility and equipment, a busy research program can be safely integrated into commercial production..

## Methods

To combine both activities in a safe way a risk-assessment should be designed to identify key control points. The use of a fully automated platform (e.g. IBA Synthera) and its disposable cassettes can prevent human errors, reduce risk of cross-contamination and improve reliability. A production schedule, training and labeling procedures should also be considered.



**IFP™ disposable cassette (Integrated Fluidic Processor)**



**Synthera® Platform designed for a GMP environment**

| RANKING SCORE | SEV = Severity of Effect on Patient  |
|---------------|--|
| 1             | Minor/No Effect - no consequence to patient safety.  |
| 2             | Mild Effect - Patient is mildly affected.  |
| 3             | Significant Effect - Severe impact on the patient. Possibility to kill or severely injure the patient. |

| Acceptable | ALARP | Not acceptable |
|------------|-------|----------------|
| 1          | 6     | 12             |
| 2          | 8     | 18             |
| 3          | 9     | 27             |
| 4          |       |                |

| RANKING SCORE | OCCURRENCE (OCC)                                  |
|---------------|---|
| 1             | Low Probability. Risk Event extremely unlikely    |
| 2             | Medium Probability-Risk Event is not frequent.    |
| 3             | High Probability. Risk Event is/could be a common |

| FMEA – Failure Mode and Effects Analysis                                   |                     |   | ANALYSIS |            |     |               |
|--|---------------------|---|----------|------------|-----|---------------|
| FAILURE MODE   | EFFECT              | CONTROLS  | SEV      | OCC (Freq) | DET | RPN S x O x D |
| Mix-ups & cross-contamination: R&D & routine production of several tracers | Business & GMP risk | Use of disposables (line clearance)<br>Appropriate training of the staff<br>Minimize human intervention<br>Well-defined planning<br>SOPs/procedures<br>Labelling ;Quality control | 3        | 2          | 1   | 6             |
| Risk of bioburden & endotoxins from starting materials                     | Business & GMP risk | Work only with qualified supplier<br>Bioburden & endotoxins control at the supplier<br>Periodic assessment of bioburden before filtration etc...                                  | 3        | 1          | 2   | 6             |

| RANKING SCORE | DETECTABILITY DESCRIPTION                    |
|---------------|--|
| 1             | Risk Event is always immediately detectable. |
| 2             | Risk Event is normally detectable.           |
| 3             | Risk Event may or may not be detected.       |

## Results



**Positronpharma** Positronpharma in Chile delivers FDG daily to several hospitals and produced > 1600 Ci FDG in the last 3 years. In addition, [<sup>18</sup>F]-FCH, NaF, FET, FLT, F-MISO, DMFP), [<sup>68</sup>Ga]-DOTA-peptides and PSMA-11 have been implemented as well as two more INDs: [<sup>18</sup>F]-PR04MZ and MHMZ were translated into clinical studies.




ICNAS a research unit of the University of Coimbra hosts a GMP-compliant PET production facility which supports clinical and pre-clinical R&D programs and supplies RPs to nearby hospitals. Currently it has five RPs authorized in the market (<sup>18</sup>F-: FDG, FCH, NaF) and <sup>68</sup>Ga-DOTA-NOC. All are produced on IBA Synthera and together represent >2000 production cycles since 2012. An extensive R&D program is in place with plans for production of other [<sup>18</sup>F]-tracers.




Synthera® set-up in PositronPharma. New Synthera® modules were purchased recently.




**Upper left:** Radiosynthesis and formulation after HPLC-purification of <sup>18</sup>F-PR04.MZ; **Center:** PET/MRI fusion, transaxial and coronal view of <sup>18</sup>F-PR04.MZ binding in striatal region, 15 min frame, 90 min p.i.; **Right:** PET/MRI fusion image of <sup>18</sup>F-MH.MZ binding to cortical and subcortical 5HT<sub>2A</sub> receptors.



<sup>68</sup>Ga-DOTA-NOC Production



R&D compounds



ICNAS production timeline

by **CYCLOTRON** VU

radiopharmaceuticals and radionuclides BV Cyclotron VU in Amsterdam is a private company with a GMP-compliant PET production facility. Since the 90's FDG is delivered for the Dutch hospitals (annual output of > 7300 patient doses (11 TBq)). <sup>18</sup>F-FCH and <sup>18</sup>F-FBB (Florbetaben) are also commercially produced with annual output >1400 (1.7 TBq) and >340 patient doses (1.1 TBq), respectively. Besides the commercial productions there is also room for an R&D program, e.g., currently improvement of the production of <sup>18</sup>F-FCH is examined.



Synthera® Synthesizer and Extension set-up at BV cyclotron for the optimization of <sup>18</sup>F-FCH process.

## Conclusion

As a conclusion, the sites described have been functioning for several years and are able to safely combine commercial daily production while keeping a highly active R&D program with more than ten different tracers developed for pre-clinical and clinical applications.