

Automated radiosynthesis of [18F]FAZA with Synthera radiochemistry box

A.G. Horti, M. Y. Kiselev, H.T. Ravert,
R. L. Wahl, R.F. Dannals


Dept. Radiology, Johns Hopkins Medicine; IBA

[18F]FAZA, 1-(5-[18F]fluoro-5-deoxy-a-D-arabinofuranosyl)-2-nitroimidazole, is an azomycin-based nucleoside for studying hypoxic tumors in patients and animals. There is a substantial interest in a simplified and automated radiosynthesis of this PET tracer. Here we present a radiosynthesis of [18F]FAZA using a new automated system called Synthera® (commercially available from Ion Beam Applications, Belgium).

All reagents and solvents used were purchased from Aldrich. The tosyl-FAZA precursor was obtained from ABX. Production of [18F]FAZA has been done via simplified method 1 using disposable integrated fluid processor (IFP™) followed by HPLC purification. The radiochemical purity and radiochemical yield of the final product were >98% and 20±4%, respectively. The time of radiosynthesis is 26 min. If required, an automatic ejection of IFP™ will allow multiple production runs per day.

Conclusion: A simple, efficient and automated radiosynthesis of [18F]FAZA has been achieved with the Synthera® (IBA) radiochemistry box.

1. Reischl, G.; Ehrlichmann, W.; Bieg, C.; Solbach, C.; Kumar, P.; Wiebe, L. I.; Machulla, H. J., *Appl Radiat Isot* 2005, 62, (6), 897-901.
2. Sorger, D.; Patt, M.; Kumar, P.; Wiebe, L. I.; Barthel, H.; Seese, A.; Dannenberg, C.; Tannappel, A.; Kluge, R.; Sabri, O. *Nucl Med Biol* 2003, 30, (3), 317-26.
3. Kumar, P.; Patt, M.; Machulla, H. J.; McEwan, A. J. B.; Wiebe, L. I., [18F]-FAZA: a putative PET radiotracer for hypoxia. *Synthesis and Applications of Isotopically Labelled Compounds, Proceedings of the International Symposium, 7th, Dresden, Germany, June 18-22, 2000*, 367-370.
4. Kumar, P.; Wiebe, L. I.; Asikoglu, M.; Tandon, M.; McEwan, A. J. *Appl Radiat Isot* 2002, 57, (5), 697-703.



Protocol

Prep HPLC: Luna C18, 50/950 ethanol/water + 0.01M NaH₂PO₄, 4 ml/min

Analytical HPLC: Supelco 50/950 ethanol/water + 0.01M NaH₂PO₄, 2 ml/min


V1=0.6 ml K222/K₂CO₃, V2=3-5 mg precursor in 0.6 ml DMSO,

V3=0.7 mL 0.2M H₃PO₄, V4=1.0 ml 0.1N NaOH,

Collect F18 from the target and start synthesis

Elute fluoride with 25 mg of K222+3 mg K₂CO₃ in 0.3mLCH₃CN+0.3mLH₂O into the reaction vessel. Dry fluoride under an argon stream and continue heating after bakeout at 120°C. Chill down to 100°C.

Add 3-5 mg of precursor in 600 µL of DMSO. Heat for 5 min at 100°C. Cool down to 25°C and add 1 mL 0.1N NaOH, Wait for 4 min at 25°C. Add 0.7mL 0.2M H₃PO₄. Transfer to the collection vial via Al₂O₃ Sep-Pak Light (washed with 10 cc H₂O). Inject the content of vial onto prep HPLC.



- Trapping of [18F]fluoride with QMA Sep-Pak
- Stripping of [18F]fluoride with K₂CO₃/Kryptofix solution
- Transfer of the reaction mixture via the Al₂O₃ Sep-Pak into the collection vial
- HPLC injection